



# Recombinant Propionate kinase (tdcD)

<b>Product Code</b>	CSB-MP841289EOD
<b>Abbreviation</b>	tdcD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8XAF1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNEFPVVLVI NCGSSSIKFS VLDASDCEVL MSGIADGINS ENAFLSVNGG EPAPLAHHSY EGALKAIAFE LEKRNLNSDV ALIGHRIAHG GSIFTESAI TDEVIDNIRR VSPLAPLHNY ANLSGIESAQ QLFPGVTQVA VFDTSFHQTM APEAYLYGLP WKYYEELGVR RYGFHGTSHR YVSLRAHSLN LAEDDSGLV VAHLGNGASI CAVRNGQSDV TSMGMTPLEG LMMGTRSGDV DFGAMSWVAS QTNQSLGDLE RVVNKESGLL GISGLSSDLR VLEKAWHEGH ERAQLAIKTF VHRIARHIAG HAASLRRLDG IIFTGGIGEN SSLIRRLVME HLAVLGVEID TEMNRSNSF GERIVSSENA RVICAVIPTN EEKMIALDAI HLGKVNAPAE FA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	tdcD
<b>Protein Names</b>	Recommended name: Propionate kinase EC= 2.7.2.15
<b>Expression Region</b>	1-402
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.