



Recombinant *Gymnogyps californianus* Serine/threonine-protein kinase mos (MOS)

Product Code	CSB-MP841925GDAQ
Abbreviation	MOS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90XV6
Product Type	Recombinant Protein
Immunogen Species	<i>Gymnogyps californianus</i> (California condor) (<i>Vultur californianus</i>)
Purity	>85% (SDS-PAGE)
Sequence	RLCLLQPLGS GGFGSVYKAT YHGATVAVKQ VKKSSKNRLA SRQSFVAELN VARLQHNNVV RVVAASTCAP ASQNSLGTII MEYVGNITLH HVIYGTGDVW RQGEDDEGGC GRKALSMEET VCYSCDIMTG LAFHLSQGIV HLDLKPANIF ITEQGVCKIG DFGCSQKLEE GLSQSPHVCQ QGGTYTHRAP ELLKGERVTA
Source	Mammalian cell
Target Names	MOS
Protein Names	Recommended name: Serine/threonine-protein kinase mos EC= 2.7.11.1 Alternative name(s): Oocyte maturation factor mos
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.