



Recombinant Danio rerio LanC-like protein 1 (lancl1)

Product Code	CSB-MP841939DIL
Abbreviation	lancl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q90ZL2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSEQRALKNP YPDYTGLGCA QDLFDMQGNL TQHFATSISS KISELLAILE NGLKNADPRD CTGYTGWAGI ALLYLHLHSV FGDPTFLQRA LDYVNRSLRS LTQRWVTFLC GDAGPLAIAA VVYHRLQKHQ ESDECLNRLL QLQPSVVQGK GRLPDELLYG RTGYLYSLIF VNQQFQQEKI PFQYIQQICD AILESGQILS QRNKIQDQSP LMYEWYQEEY VGAAHGLSGI YYYLMQPGLV AGQDRVFSLV KPSVNYVCQL KFPSGNYAPC VGDARDLLVH WCHGSPGVIY MLIQAFKVFV VRQYLEDALQ CGEVIWQRGL LKKGYGLCHG AAGNAYGFLA LYKITQDPKH LYRACMFADW CMNYGRHGCR TPDTPFSLFE GMAGTIYFLA DLLQPARAKF PCFEV
Source	Mammalian cell
Target Names	lancl1
Protein Names	Recommended name: LanC-like protein 1
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.