



Recombinant Rat Cryptochrome-2 (Cry2)

Product Code	CSB-MP842112RA
Abbreviation	Cry2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q923I8
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAVVAAT VPVQSMGADG ASSVHWFRKG LRLHDNPALL AAVRGARCVR CVYILDPWFA ASSSVGINRW RFLQSLLEDL DTSRKLNSR LFVVRGQPAD VFPRLFKEWG VTRLTFEYDS EPFGKERDAA IMKMAKEAGV EVTTENSHTL YDLDRIIELN GQKPPLTYKR FQALISRMEL PKKPVGAVSS QHMCNCRAEI QENHDDTYGV PSLEELGFPT EGLGPAVWQG GETEALVRLD KHLERKAWVA NYERPRMNAN SLLASPTGLS PYLRFGLCLSC RLFYYRLWDL YRKVKRNSTP PLSLFGQLLW REFFYTAATN NPRFDRMEGN PICIQIPWDR NPEALAKWAE GKTGFPWIDA IMTQLRQEGW IHHLARHAVA CFLTRGDLWV SWESGVRVFD ELLLDADFSV NAGSWMWLSC SAFFQQFFHC YCPVGFGRRT DPSGDYIRRY LPKLGFPSPR YIYEPWNAPE SVQKAANCII GVDYPRPIVN HAETSRLNIE RMKQIQQLS RYRGLCLWAS VPSCVEDLSH PVAEPGSSQA GSISNTGPRP LSSGPASPKR KLEAAEPPG EELSKRARVT VTQMPAQEPP SKDS
Source	Mammalian cell
Target Names	Cry2
Protein Names	Recommended name: Cryptochrome-2
Expression Region	1-594
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.