



# Recombinant Human Acid sphingomyelinase-like phosphodiesterase 3a (SMPDL3A)

<b>Product Code</b>	CSB-BP842119HU
<b>Abbreviation</b>	SMPDL3A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92484
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LPVAPAGG RNPPPAIGQF WHVTDLHLDP TYHITDDHTK VCASSKGANASNPDPFGDVL CDSPLYQLILS AFDFIKNSGQ EASFMIWTGD SPPHVPVPELSTDTVINVIT NMTTTIQSLF PNLQVFPALG NHDYWPQDQL PVRTSKVYNAVANLWKPWLD EEAISTLRKG GFYSQKVTTN PNLRIISLNT NLYYGPNIMTLNKTDPANQF EWLESTLNNS QQNKEKVYII AHVPVGYLPS SQNITAMREY YNEKLIDIFQ KYSDVIAGQF YGHTRDSIM VLSDKKGSPV NSLFVAPAVT PVKSVLEKQT NNPGIRLFQY DPRDYKLLDM LQYYLNLTEA NLKGESIWKLEYILTQTYDI EDLQPESLYG LAKQFTILDS KQFIKYNYF FVSYDSSVTC DKTCKAFQIC AIMNLDNISY ADCLKQLYIK HNY
<b>Source</b>	Baculovirus
<b>Target Names</b>	SMPDL3A
<b>Protein Names</b>	Recommended name: Acid sphingomyelinase-like phosphodiesterase 3a Short name= ASM-like phosphodiesterase 3a EC= 3.1.4.-
<b>Expression Region</b>	23-453
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.