



# Recombinant Rat Endothelial transcription factor GATA-2 (Gata2)

<b>Product Code</b>	CSB-MP842125RA
<b>Abbreviation</b>	Gata2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q924Y4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVAPEQPRW MAHPAVLNAQ HPDSHHPGLA HNYMEPAQLL PPDEVDFVFN HLDSQGNPYY ANPAHARARV SYSPAHARLT GGQMCRPHLL HSPGLPULDG GKAALSAAAA HHHSPWTVSP FSKTPLHPSA AGAPGGPLSV YPGAAGGSGG GSGSSVASLT PTAHSGSHL FGFPPTPPKE VSPDPSTTGA ASPASPSAGG SVARGEDKDG VKYQVSLSES MKMEGGSPLR PGLAPMGTQP ATHHPIPTY SYVPASAHEY GSGLFHPPGGF LGGPASSFTP KQRSKARSCS EGRECVNCGA TATPLWRRDG TGHYLCNACG LYHKMNGQNR PLIKPKRRLS AARRAGTCCA NCQTTTTTLW RRNANGDPVC NACGLYYKLH NVNRPLTMKK EGIQTRNRKM SSKSKKSKKG AECFEELSKC MQEKSSPFSA AALAGHMAPV GHLPPFSHSG HILPTPTPIH PSSSLFSGHP HPSSMVTAMG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Gata2
<b>Protein Names</b>	Recommended name: Endothelial transcription factor GATA-2 Alternative name(s): GATA-binding protein 2
<b>Expression Region</b>	1-480
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.