



# Recombinant Human UBX domain-containing protein 4 (UBXN4)

<b>Product Code</b>	CSB-EP842131HU
<b>Abbreviation</b>	UBXN4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92575
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLWFQGAIPA AIATAKRSGA VFVVFVAGDD EQSTQMAASW EDDKVTEASS NSFVAIKIDT KSEACLQFSQ IYPVVCVPSS FFIGDSGIPL EVIAGSVSAD ELVTRIHKVR QMHLLKSETS VANGSQSESS VSTPSASFEP NNTCENSQSR NAELCEIPT SDTKSDTATG GESAGHATSS QEPSGCSDQR PAEDLNIRVE RLTKKLEERR EEKRKEEQEQR EIKKEIERRK TGKEMLDYKR KQEEELTKRM LEERNREKAE DRAARERIKQ QIALDRAERA ARFAKTKEEV EAKAAALLA KQAEMEVKRE SYARERSTVA RIQFRLPDGS SFTNQFPSDA PLEEARQFAA QTVGNTYGNF SLATMFPRRE FTKEDYKKKL LDLELAPSAS VLLPAGRPT ASIVHSSSGD IWTLTGTVLY PFLAIWRLIS NFLFSNPPPT QTSVRVTSSE PPNPASSSKS EKREPVRKRV LEKRGDDFKK EGKIYRLRTQ DDGEDENNTW NGNSTQQM
<b>Source</b>	E.coli
<b>Target Names</b>	UBXN4
<b>Protein Names</b>	Recommended name: UBX domain-containing protein 4 Alternative name(s): Erasin UBX domain-containing protein 2
<b>Expression Region</b>	1-508
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.