



Recombinant Human UBX domain-containing protein 4 (UBXN4)

Product Code	CSB-MP842131HU
Abbreviation	UBXN4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92575
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLWFQGAIPA AIATAKRSGA VFVVFVAGDD EQSTQMAASW EDDKVTEASS NSFVAIKIDT KSEACLQFSQ IYPVVCVPSS FFIGDSGIPL EVIAGSVSAD ELVTRIHKVR QMHLLKSETS VANGSQSESS VSTPSASFEP NNTCENSQSR NAELCEIPT SDTKSDTATG GESAGHATSS QEPGCSQDQR PAEDLNIRVE RLTKKLEERR EEKRKEEQEQR EIKKEIERRK TGKEMLDYKR KQEEELTKRM LEERNREKAE DRAARERIKQ QIALDRAERA ARFAKTKEEV EAKAAALLA KQAEMEVKRE SYARERSTVA RIQFRLPDGS SFTNQFPSDA PLEEARQFAA QTVGNTYGNF SLATMFPRRE FTKEDYKKKL LDLELAPSAS VLLPAGRPT ASIVHSSSGD IWTLTGTVLY PFLAIWRLIS NFLFSNPPPT QTSVRVTSSE PPNPASSSKS EKREPVRKRV LEKRGDDFKK EGKIYRLRTQ DDGEDENNTW NGNSTQQM
Source	Mammalian cell
Target Names	UBXN4
Protein Names	Recommended name: UBX domain-containing protein 4 Alternative name(s): Erasin UBX domain-containing protein 2
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.