



Recombinant Methylthioribose-1-phosphate isomerase (C01G10.9)

Product Code	CSB-YP842357CXY
Abbreviation	C01G10.9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93169
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTSKINKIVN GPIPDSFLFD ESKRLDSLKF DGTNLEVLDQ LLLPHEFKYI PVEGVSDAFA VIKSMQVRGA PLIAVVGSLG LLEIQKASE LDESIIQKI NFLISSRPTA VDLRNSLNGL KPILESQDYS DVVKLEKCRS YLLNVYTDEK LQNRILVWNA YQELLSAFPD KEKLTVM TIC NTGSLATISW GTALGVIRAL HSENRLKLVY VLETRPYNQG IRLTSIELLH GEVPFKLITD SMAAWAMKNH QVDCVLTGAD NVARNGDTAN KIGTYMLAVL CKHHNINFYP VVPFTTINKN IATGEEIKIE ERPSAEMLRV NGVLIGNSEC PVWNPAPFDVT PAHLITKILT DFGNWPPEML EQQIPK
Source	Yeast
Target Names	C01G10.9
Protein Names	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S-methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.