



Recombinant Arabidopsis thaliana Phosphoenolpyruvate carboxylase kinase 2 (PPCK2)

Product Code	CSB-EP842406DOA
Abbreviation	PPCK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93VK0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTREFELENN YQLCDEIGRG RFGTITRCFS PATKEYACK TIDKRVLIDA LDRECIETEP RIMAMLPPHP NIIRIFDLYE TEDSLAIVME LVDPPMTIYD RLISAGGRLS ESESASYAKQ ILSALAHCHR CDVVHRDVKP DNVLVDLVSG GVKLCDFGSA VWLGGETAEG VVGTPYYVAP EVVMGRKYDE KVDIWSAGVV IYTMLAGEPP FNGETAEDIF ESILRGNLRF PPKKFGSVSS EAKDLLRKMI CRDVSRRFSA EDALRHSWMM NVGNLQSN
Source	E.coli
Target Names	PPCK2
Protein Names	Recommended name: Phosphoenolpyruvate carboxylase kinase 2 Short name= AtPPCK2 EC= 2.7.11.1
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.