



Recombinant Arabidopsis thaliana Adenylate isopentenyltransferase 3, chloroplastic (IPT3)

Product Code	CSB-YP842409DOA
Abbreviation	IPT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93WC9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	RLSVD IATRFRAEII NSDKIQVHQG LDIVTNKITS EESCGVPHHL LGVLPEADL TAANYCHMAN LSIESVLNRG KLPIIVGGSN SYVEALVDDK ENKFRSRYDC CFLWVDVALP VLHGFVSERV DKMVESGMVE EVREFFDFSN SDYSRGIKKA IGFPEFDRFF RNEQFLNVED REELLSKYLE EIKRNTFELA CRQREKIERL RKVKKWSIQR VDATPVFTKR RSKMDANVAW ERLVAGPSTD TVSRFLLDIA SRRPLVEAST AVAAAMEREL SRC
Source	Yeast
Target Names	IPT3
Protein Names	Recommended name: Adenylate isopentenyltransferase 3, chloroplastic Short name= AtIPT3 EC= 2.5.1.n4 Alternative name(s): Adenylate dimethylallyltransferase 3 Cytokinin synthase 3
Expression Region	56-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.