



Recombinant Arabidopsis thaliana Probable beta-1,4-xylosyltransferase IRX10L (IRX10L), partial

Product Code	CSB-BP842430DOA
Abbreviation	IRX10L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q940Q8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Source	Baculovirus
Target Names	IRX10L
Protein Names	Recommended name: Probable beta-1,4-xylosyltransferase IRX10L EC= 2.4.2.- Alternative name(s): Glucuronoxylan glucuronosyltransferase 1 Short name= AtGUT1 Glucuronoxylan glucuronosyltransferase 2 Short name= AtGUT2 Pro
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.