



Recombinant Arabidopsis thaliana F-box protein At5g03970 (At5g03970)

Product Code	CSB-YP842463DOA
Abbreviation	At5g03970
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94A64
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSCSARLKTN LGMTSALSTH EVLNSNDTMC EILILLPPET IYKLILVSKR WLEIIASPCF RHTYLAKWKP NFELIGFFVC SSMYLGRRID GTRRPRAEPS LPLLSTSSIG DEIESSGILK KLGYYIDSAD GLLLCGRHPK AYYLWDPSTQ KQQQLPRPRV HFEELCMSLI SEDSPDQGFS YKVVRAECVA FPNSTKVKV ETFSSKTSTW SYSELICLEP LSLSPWTPGR VIKGVVYWHA RGGKIAIYDS NSEEKRIDVI KLPKTYDYDE QVLGETDDGS LQYGWSNKS V MEVWKLEKVG NVLEWNLLFK VNFKAMWKMN QAVATRFSTW TKETQLLALF NQKSDLVYIR CDTQIFIYHT ETKRVEVVQY QGRKSTLVWD FNKVVYPYFRR TWPCSPLEFEN NAMMIDTPQC NASSAGSN
Source	Yeast
Target Names	At5g03970
Protein Names	Recommended name: F-box protein At5g03970
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.