



Recombinant Arabidopsis thaliana Glutathione S-transferase F8, chloroplastic (GSTF8)

Product Code	CSB-EP842594DOA-B
Abbreviation	GSTF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96266
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	A SIKVHGVPMS TATMRVLATL YEKDLQFELI PVDMRAGAHK QEAHLALNPF GQIPALEGD LTLFESRAIT QYLAEYSEK GEKLISQDCK KVKATTNVWL QVEGQQFDPN ASKLAFERVF KGMFGMTTDP AAVQELEGKL QKVLVDVYEAR LAKSEFLAGD SFTLADLHHL PAIHYLLGTD SKVLFDSRPK VSEWIKKISA RPAWAKVIDL QKQ
Source	E.coli
Target Names	GSTF8
Protein Names	Recommended name: Glutathione S-transferase F8, chloroplastic Short name= AtGSTF8 EC= 2.5.1.18 Alternative name(s): AtGSTF5 GST class-phi member 8 Glutathione S-transferase 6
Expression Region	50-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.