



Recombinant Mediator of RNA polymerase II transcription subunit 18 (mdt-18)

Product Code	CSB-YP842617CXY
Abbreviation	mdt-18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q966M5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MDNLVDIDED GTQEKLSVNS TPYQTQECVL YGSIFAKNVP DLERRLAGLC DPGCEEFFEH EMSFSLRSSG IMDIKTDIKL RRRFRTEHRI QNYWQLKYIG VPEPDQKCPT IVRKEISLV HSQDMMTYAK SLGLRMDY EY MTQGKLWTKG NIKILHSILT STSRAGTYDS NSIKSVSDSA LVEISISLPE SAEYMPAAKS LRDFADQLMP LVNMEKVDYW KKMFPQAAP RR
Source	Yeast
Target Names	mdt-18
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 18 Alternative name(s): Mediator complex subunit 18
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.