



# Recombinant Human Transcription elongation factor A protein-like 3 (TCEAL3)

<b>Product Code</b>	CSB-EP842623HU-B
<b>Abbreviation</b>	TCEAL3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q969E4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEKPYNKNEG NLENEGKPED EVEPDDEGKS DEEEKPDVEG KTECEGKRED EGEPGDEGQL EDEGSQEKQG RSEGEGKPQG EGKPASQAKP ESQPRAAEKR PAEDYVPRKA KRKTDRGTDD SPKDSQEDLQ ERHLSSEEMM RECGDVSRAQ EELRKKQKMG GFHWMQRDVQ DPFAPRGQRG VRGVRGGGRG QRGLHDIPYL
<b>Source</b>	E.coli
<b>Target Names</b>	TCEAL3
<b>Protein Names</b>	Recommended name: Transcription elongation factor A protein-like 3 Short name= TCEA-like protein 3 Alternative name(s): Transcription elongation factor S-II protein-like 3
<b>Expression Region</b>	1-200
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .