



# Recombinant Human SH3 domain-containing kinase-binding protein 1 (SH3KBP1)

<b>Product Code</b>	CSB-EP842640HU-B
<b>Abbreviation</b>	SH3KBP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96B97
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVEAIVEFDY QAQHDELTI SVGEIITNIR KEDGGWEGQ INGRRLFPD NFVREIKKEM KKDPLTNKAP EKPLHEVPSG NSLLSSETIL RTNKRGERRR RRCQVAFSYL PQNDDELELK VGDIIEVVEE VEEGWEGVL NGKTGMFPSN FIKELSGESD ELGISQDEQL SKSSLRETTG SESDGGDSSS TKSEGANSTV ATAAIQPKKV KGVGFGDIFK DKPIKLRPRS IEVENDFLPV EKTIGKKLPA TTATPDSSKT EMDSRTKSKD YCKVIFPYEA QNDELTIKE GDIVTLINKD CIDVGWWEGE LNRRRGVFPD NFVKLLPPDF EKEGNRPKKP PPSAPVIKQ GAGTTERKHE IKKIPPERPE MLPNRTEEKE RPEREPKLDL QKPSVPAIPP KKPRPPKTNS LSRPGALPPR RPERVVGPLT HTRGDSPKID LAGSSLSGIL DKDLSDRSND IDLEGFDSVV SSTEKLSHPT TSRPKATGRR PPSQSLTSSS LSSPDIFDSP SPEEDKEEHI SLAHRGVDAS KTKSTVTIS QVSDNKASLP PKPGTMAAGG GGPAPLSSAA PSPLSSSLGT AGHRANSPSL FGTEGKPKME PAASSQAAVE ELRTQVREL RSIETMKDQQ KREIKQLLSE LDEEKIRLR LQMEVNDIKK ALQSK
<b>Source</b>	E.coli
<b>Target Names</b>	SH3KBP1
<b>Protein Names</b>	Recommended name: SH3 domain-containing kinase-binding protein 1 Alternative name(s): CD2-binding protein 3 Short name= CD2BP3 Cbl- interacting protein of 85 kDa Human Src family kinase-binding protein 1 Short name= HSB-1
<b>Expression Region</b>	1-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.