



Recombinant Human Protein Spindly (CCDC99)

Product Code	CSB-EP842662HU-B
Abbreviation	CCDC99
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EA4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEADIITNLR CRLKEAEEER LKAAQYGLQL VESQNELQNNQ LDKCRNEMMT MTESYEQEKY TLQREVELKS RMLESLSCEC EAIKQQQKM HLEKLEEQLSR SHGQEVNELK TKIEKLKVEL DEARLSEKQL KHQVDHQKEL LSCKSEELRV MSERVQESMS SEMALQIEL TEMESMKTTL KEEVNELQYR QEQLLELLITN LMRQVDRLEKE EKEEREKEAV SYYNALAEKAR VANQDLQVQL DQALQQALDP NSKGNLSLFAE VEDRRAAMER QLISMVKYQ SLKKQNVFNR EQMQRMKLQI ATLLQMKGSQ TEFEQQERLL AMLEQKNGEI KHLLGEIRNL EKFKNLYDSM ESKPSVDSGT LEDNTYYTDL LQMKLDNLNK EISTKGELS IQRMKALFES QRALDIERKL FANERCLQLS ESENMKLRAK LDELKLKYEP EETVEVPVLK KRREVLVVDI TTAKDACVNN SALGGEVYRL PPQKEETQSC PNSLEDNNLQ LEKSVSIYTP VVSLSPHKNL PVDMLKKEK KCVKLIGVPA DAEALSERSG NTPNSPRLAA ESKLQTEVKE GKETSSKLEK ETCKKLHPIL YVSSKSTPET QCPQQ
Source	E.coli
Target Names	SPDL1
Protein Names	Recommended name: Protein Spindly Short name= hSpindly Alternative name(s): Arsenite-related gene 1 protein Coiled-coil domain-containing protein 99 Rhabdomyosarcoma antigen MU-RMS-40.4A
Expression Region	1-605
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.