



Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 (ZHX1-C8orf76)

Product Code	CSB-EP842663HU-B
Abbreviation	ZHX1-C8orf76
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96EF9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLRKLWQWFY EETESSDDVE VLTLKKFKGD LAYRRQEYQK ALQEYSSISE KLSSTNFAMK RDVQEGQARC LAHLGRHMEA LEIAANLENK ATNTDHLTTV LYLQLAICSS LQNLEKTIFC LQKLISLHPF NPWNWGKLAE AYLNLGPALS AALASSQKQH SFTSSDKTIK SFFPHSGKDC LLCFPETLPE SSLFSVEANS SNSQKNEKAL TNIQNCMAEK RETVLIETQL KACASFIRTR LLLQFTQPQQ TSFALERNLR TQQEIEDKMK GFSFKEDTLL LIAEVSVSLG FM
Source	E.coli
Target Names	ZHX1-C8orf76
Protein Names	Recommended name: Zinc fingers and homeoboxes protein 1, isoform 2 Alternative name(s): ZHX1-C8orf76 readthrough transcript protein
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.