



# Recombinant Human ELMO domain-containing protein 3 (ELMOD3)

<b>Product Code</b>	CSB-EP842670HU-B
<b>Abbreviation</b>	ELMOD3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96FG2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNEKSCSFHS KEELRDGQGE RLSAGYSPSY DKDKSVLAFR GIPISLKNH GILQALTTEA YEWEPVVST EVVRAQEEWE AVDTIQPETG SQASSEQPGQ LISFSEALQH FQTVDLSPFK KRIQPTIRRT GLAALRHLYF GPPKLHQRLR EERDLVLTIA QCGLDSQDPV HGRVLQTIYK KLTGSKFDCA LHGNHWEDLG FQGANPATDL RGAGFLALLH LLYLVMSKT LPMAQEIFRL SRHHIQQFPF CLMSVNITHI AIQALREECL SRECNRQQKV IPVVNSFYAA TFLHLAHVWR TQRKTISDSG FVLKELEVLA KKSPPRLLKT LELYLARVSK GQASLLGAQK CYGPEAPPFK DLTFTGESDL QSHSSEGVWL I
<b>Source</b>	E.coli
<b>Target Names</b>	ELMOD3
<b>Protein Names</b>	Recommended name: ELMO domain-containing protein 3 Alternative name(s): RNA-binding motif and ELMO domain-containing protein 1 RNA-binding motif protein 29 RNA-binding protein 29
<b>Expression Region</b>	1-381
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.