



Recombinant Human Zinc finger C2HC domain-containing protein 1A (ZC2HC1A)

Product Code	CSB-EP842680HU-B
Abbreviation	ZC2HC1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96GY0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGLEENGGV VQVGELLPCK ICGRTFFPVA LKKHGPICQK TATKKRKTFD SSRQRAEGTD IPTVKPLKPR PEPPKKPSNW RRKHEEFIAT IRAAKGLDQA LKEGGKLPPP PPSYDPDYI QCPYCQRRFN ENAADRHINF CKEQAARISN KGFSTDTKG KPTSRTQVYK PPALKKSNP GTASSGSSRL PQPSGAGKTV VGVPSGKVSS SSSSLGNKLQ TLSPSHKGIA APHAGANVKP RNSTPPSLAR NPAPGVLTKN RKTYTESYIA RPDGDCASSL NGGNIKGIEG HSPGNLPKFC HECGTKYPVE WAKFCCECGI RRMIL
Source	E.coli
Target Names	ZC2HC1A
Protein Names	Recommended name: Zinc finger C2HC domain-containing protein 1A
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.