



Recombinant Human HRAS-like suppressor 5 (HRASLS5)

Product Code	CSB-MP842699HU
Abbreviation	HRASLS5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96KN8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGLSPGAEGE YALRLPRIPP PLPKPASRTA STGPKDQPPA LRRSAVPHSG LNSISPLELE ESVGFAALVQ LPAKQPPPGT LEQGRSIQQG EKAVVSLETT PSQKADWSSI PKPENEGKLI KQAAEGKPRP RPDGLIEIFR IGYEHWAIYV EDDCVVHLAP PSEEFVSGSI TSIFSNRAVV KYRLEDVLH GCSWKVNNKL DGTYLPLPVD KIIQRTKKMV NKIVQYSLIE GNCEHFVNLG RYGVPRSQQV EHALMEGAKA AGAVISAVVD SIKPKPITA
Source	Mammalian cell
Target Names	HRASLS5
Protein Names	Recommended name: HRAS-like suppressor 5 Alternative name(s): H-rev107-like protein 5
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.