



Recombinant Human PWWP domain-containing protein 2A (PWWP2A)

Product Code	CSB-MP842723HU
Abbreviation	PWWP2A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96N64
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAVAAEAAA TAASPGEGGA GEAPEMEPI PGSEAGTDPL PVTATEASVP DGETDGQQSA PQADEPPLPP PPPPGELAR SPEAVGPELE AEEKLSVRVA ESAAAAPQGG PELPPSPASP PEQPPAPEER EEPPLPQPVA PALVPPAGGD STVSQIPGS EVRVTLDHII EDALVVSFRF GEKLFSGVLM DLSKRFGPHG IPVTVPKRE YKDKPEAMPL QSNTFQEGTE VKCEANGAVP DDPSPVPHPE LSLAESLWTS KPPPLFHEGA PYPPPLFIRD TYNQSIPQPP PRKIKRPKRK MYREEPTSIM NAIKLRPRQV LCDKCKNSVV AEKKEIRKGS SATDSSKYED KKRRNESVTT VNKKLKT DHK VDGKNQNESQ KRNAVVKVSN IAHSRGRVVK VSAQANTSKA QLSTKKVLQS KNMDHAKARE VLKIAKEKAQ KKQNETSTSK NAHSKVHFTR RYQNPSSGSL PPRVRLKPQR YRNEENDSSL KTGLEKMRSR KMAPKPQSRC TSTRSAGEAP SENQSPSKGP EEASSEVQDT NEVHVP GDQD EPQTLGKKGS KNNISVYMTL NQKKSDSSSA SVCSIDSTDD LKSSNSECSS SESFDFPPGS MHAPSTSSTS SSSKEEKKLS NSLKMKVFSK NVSKCVTPDG RTICVGDIW AKIYGFPWWP ARILTITVSR KDNGLLVRQE ARISWFGSPT TSFLALSQLS PFLENFQSRF NKKRKGLYRK AITEAAKAAK QLTPEVRALL TQFET
Source	Mammalian cell
Target Names	PWWP2A
Protein Names	Recommended name: PWWP domain-containing protein 2A
Expression Region	1-755
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.