



Recombinant Human Zinc finger imprinted 3 (ZIM3)

Product Code	CSB-YP842734HU
Abbreviation	ZIM3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96PE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNNSQGRVTF EDVTVNFTQG EWQRLNPEQR NLYRDVMLEN YSNLVSVGQG ETTKPDVILR LEQGKEPWLE EEEVLGSGRA EKNGDIGGQI WKPKDVKESL AREVPSINKE TLTTQKGVEC DGSKKILPLG IDDVSSLQHY VQNNSHDDNG YRKLVGNNPS KFGVQQLKCN ACRKLFSSKS RLQSHLRRHA CQKPFEC HSC GRAFGEKWKL DKHQKTHAEE RPYKCENCGN AYKQKSNLFQ HQKMHTKEKP YQCKTCGKAF SWKSSCINHE KIHNAKKS YQ CNECEKSFRQ NSTLIQHKKV HTGQKPFQCT DCGKAFIYKS DLVKHQRIHT GEKPYKCSIC EKAFSQKSNV IDHEKIHTGK RAYECDLCGN TFIQKKNLIQ HKKIHTGEKP YECNRCGKAF FQKSNLHSHQ KTHSGERTYR CSECGKTFIR KLNL SLHKKT HTGQKPYGCS ECGKAFADRS YLVRHQKRIH SR
Source	Yeast
Target Names	ZIM3
Protein Names	Recommended name: Zinc finger imprinted 3 Alternative name(s): Zinc finger protein 657
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.