



Recombinant Human UAP56-interacting factor (FYTTD1)

Product Code	CSB-MP842744HU
Abbreviation	FYTTD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QD9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNRFGTRLVG ATATSSPPPK ARSNENLDKI DMSLDDIIKL NRKEGKKQNF PRLNRRLLQQ SGAQQFRMRV RWGIQQNSGF GKTSLNRRGR VMPGKRRPNG VITGLAARKT TGIRKGISPM NRPPLSDKNI EQYFPVLKRK ANLLRQNEGQ RKPVAVLKRP SQLSRKNNIP ANFTRSGNKL NHQKDTRQAT FLFRRGLKVQ AQLNTEQLLD DVVAKRTRQW RTSTTNGGIL TVSIDNPGAV QCPVTQKPRL TRTAVPSFLT KREQSDVKKV PKGVPLQFDI NSVGKQTGMT LNERFGILKE QRATLTYNKG GSRFVTVG
Source	Mammalian cell
Target Names	FYTTD1
Protein Names	Recommended name: UAP56-interacting factor Alternative name(s): Forty-two-three domain-containing protein 1 Short name= Protein 40-2-3
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.