



# Recombinant *Sulfolobus solfataricus* Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-BP842962FPM
<b>Abbreviation</b>	srp54
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q97ZE7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Sulfolobus solfataricus</i> (strain ATCC 35092 / DSM 1617 / JCM 11322 / P2)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLENIRDAVR KFLTGSTPYE KAVDEFIKDL QKSLISSDVN VKLVFSLTAK IKERLNKEKP PSVLERKEWF ISIVYDELSK LFGGDKEPNV NPTKLPFIIM LVGVQGSQKT TTAGKLAYFY KKRQYKVGVLV AADVYRPAAY DQLLQLGNQI GVQVYGEPPN QNPIEIAKKG VDIFVKNKMD IIIVDTAGRH GYGEETKLE EMKEMYDVLK PDDVILVIDA SIGQKAYDLA SRFHQASPIG SVIITKMDGT AKGGGALSAV VATGATIKFI GTGEKIDELE TFNAKRFVSR ILGMGDIESI LEKVKGLEEY DKIQQKMEDV MEGKGKLTLR DVYAQIIALR KMGPLSKVLQ HIPGLGIMLP TPSEDQLKIG EEKIRRWLAA LNSMTYKELE NPNIIDKSRM RRIAEGSGLE VEEVRELLEW YNNMNRLLKM VKRRRRGNIDK LFGGKIG
<b>Source</b>	Baculovirus
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.