



Recombinant Danio rerio Homeobox protein Dlx4a (dlx4a)

Product Code	CSB-MP842995DIL
Abbreviation	dlx4a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q98879
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTMTSLSESL EASDPSKSAF LEFSHGYQSH QQHSPGVSHA HYPVHGLHQG AHSQYDAAFS PGAASYSRPL AYHYSTAHHH PGAYLPYQHN SAVGYSRVED ADSEKQSSIE SGEIRLNGKG KKIRKPRTIY SSLQLQALNQ RFQQTQYLAL PERADLAAKL GLTQTQVKIW FQNKRSKYKK IMKHGSSGPE GEHLQAASAS GAPCSPGMPP LWDVSMPSKG APIHSGGYMN SFGHWYSGHH QDPMARTQMM
Source	Mammalian cell
Target Names	dlx4a
Protein Names	Recommended name: Homeobox protein Dlx4a Alternative name(s): DLX-8 Distal-less homeobox protein 4a
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.