



Recombinant *Acanthocheilonema viteae* Fatty-acid and retinol-binding protein 1 (far-1)

Product Code	CSB-BP843128AVT
Abbreviation	far-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8MZJ8
Product Type	Recombinant Protein
Immunogen Species	<i>Acanthocheilonema viteae</i> (Filarial nematode worm) (<i>Dipetalonema viteae</i>)
Purity	>85% (SDS-PAGE)
Sequence	NVVP FSMSNIPEEY KEFIPEEVKN FYKNLTQEDR QILRELASKH ATFTNEDAAL EALKNKSDKL YQKAVELRNF VKAKIDSLKP DAKAFVDEII AKVRSRLED GQKLDVEKLN QAARDIIAKY EALNEETKEE LKAPFPNTTK IITNEKFTRI ANSFLQKN
Source	Baculovirus
Target Names	far-1
Protein Names	Recommended name: Fatty-acid and retinol-binding protein 1 Alternative name(s): Av-FAR-1 Av20
Expression Region	17-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.