



Recombinant Human Protein canopy homolog 4 (CNPY4)

Product Code	CSB-EP843134HU
Abbreviation	CNPY4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N129
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GMLKEEDDD TERLPSKCEV CKLLSTELQA ELSRTGRSRE VLELGQVLDT GKRKRHVPYS VSETRLEEAL ENLCERILDY SVHAERKGS L RYAKGQSQTM ATLKGLVQKG VKVDLGIPLE LWDEPSVEVT YLKKQCETML EEFEDIVGDW YFHHQEQLPQ NFLCEGHVLP AAETAQLQET WTGKEITDGE EKTEGEEEEQE EEEEEEEEEG GDKMTKTGSH PKLDREDL
Source	E.coli
Target Names	CNPY4
Protein Names	Recommended name: Protein canopy homolog 4
Expression Region	22-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.