



Recombinant Human Protein FAM63A (FAM63A)

Product Code	CSB-EP843272HU
Abbreviation	FAM63A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N5J2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEYHQPEDPA PGKAGTAEAV IPENHEVLAG PDEHPQDTDA RDADGEARER EPADQALLPS QCGDNLESPL PEASSAPPGP TLGTLPEVET IRACSMQPQL PQSPRTRQPE PDFYCVKWIP WKGEQTPIIT QSTNGPCPLL AIMNILFLQW KVKLPPQKEV ITSDELM AHL GNCLLSIKPQ EKSEGLQLNF QQNVDDAMTV LPKLATGLDV NVRFTGVSD F EYTPECSVFD LLGIPLYHGW LVDPQSPEAV RAVGKLSYNQ LVERIITCKH SSDTNLVTEG LIAEQFLETT AAQLTYHGLC ELTAAAKEGE LSVFFRNNHF STMTKHKSHL YLLVTDQGFL QEEQVWVESL HNVDGDSCFC DSDFHLSHSL GKGPAGAE GGS GSPETQLQVD QDYLIALLSLQ QQQPRGPLGL TDLELAQQLQ QEEYQQQQA QPVRMRTRVL SLQGRGATSG RPAGERRQRP KHESDCILL
Source	E.coli
Target Names	MINDY1
Protein Names	Recommended name: Protein FAM63A
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.