



Recombinant Human TBC1 domain family member 19 (TBC1D19)

Product Code	CSB-EP843274HU-B
Abbreviation	TBC1D19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N5T2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLQEESDLSL IIAQIVQKLK GSNLYSQLER QAWASLQRPE IKLESLKEDI KEFFKISGWE KKLQNAVYSE LSVFPLPSHP AAPPEHLKEP LVYMRKAQGS WEKRILKSLN SMCTELSIPL ARKRPVGEQK ELLNKWNEMG TDEPDLSLFR PVYAPKDFLE VLINLRNPNY ENGDSLSEFRT HLGLIQVPLK VKDIPELKEC FVELGLNIGQ LGIDDSTQVP PELFENEHVR IGQKVLAEQD SAAAQQYIRQ GSPTALRAEL WALILNISSQ PEDVLYEQL KTNVIQHDLL VDSLIIKDVK LTASNDDYYF VFEDYLYQVL LCFSRDTSVL SHFAFNSASP PKSYIRGKLG LEEYAVFYPP NGVIPFHGFS MYVAPLCFLY HEPKLYQIF REMYVRRFFR LHSISSHPSG IVSLCLLFET LLQTYLPQLF YHLREIGAQP LRISFKWMVR AFSGYLATDQ LLLLWDRILG YNSLEILAVL AAVFAFRAV NLMEVTSALAA AEAVLADLST LKVMPLQLIF LFATVT
Source	E.coli
Target Names	TBC1D19
Protein Names	Recommended name: TBC1 domain family member 19
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.