



Recombinant Human RNA polymerase II-associated factor 1 homolog (PAF1)

| | |
|--------------------------|--|
| Product Code | CSB-EP843285HU |
| Abbreviation | PAF1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8N7H5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | <p>MAPTIQTQAQ REDGHRPNSH RTLPERSGVV CRVKYCNSLP DIPFDPKFIT YFPDQNRVQ YKATSLEKQH KHDLLTEPDL GVTIDLINPD TYRIDPNVLL DPADEKLLLE EIQAPTSSKR SQQHAKVVPW MRKTEYISTE FNRYGISNEK PEVKIGVSVK QQFTEEEIYK DRDSQITAIE KTFEDAQKSI SQHYSKPRVT PVEVMPVFPD FKMWINPCAQ VIFDSDPAPK DTSGAAALEM MSQAMIRGMM DEEGNQFVAY FLPVEETLKK RKRQDQEEEMD YAPDDVYDYK IAREYNWNVK NKASKGYEEN YFFIFREGDG VYYNELETRV RLSKRRAKAG VQSGTNALLV VKHRDMNEKE LEAQEARKAQ LENHEPEEEEE EEMETEEKE AGGSDEEQEK GSSSEKEGSE DEHSGSESER EEGDRDEASD KSGSGEDESS EDEARAARDK EEIFGSDADS EDDADSDDDED RGQAQGGSDN DSDSGSNGGG QRSRSHRSA SPFPSGSEHS AQEDGSEAAA SDSSEADSDS D</p> |
| Source | E.coli |
| Target Names | PAF1 |
| Protein Names | Recommended name: RNA polymerase II-associated factor 1 homolog Short name= hPAF1 Alternative name(s): Pancreatic differentiation protein 2 |
| Expression Region | 1-531 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.