



Recombinant Human RNA polymerase II-associated factor 1 homolog (PAF1)

Product Code	CSB-MP843285HU
Abbreviation	PAF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N7H5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAPTIQTQAQ REDGHRPN SH RTLPERSGVV CRVKYCNSLP DIPFDPKFIT YFPDQNRVQ YKATSLEKQH KHDLLTEPDL GVTIDLINPD TYRIDPNVLL DPADEKLL EE EIQAPTS SKR SQQHAKVVPW MRKTEYISTE FNRYGISNEK PEVKIGVSVK QQFTEEEIYK DRDSQITAIE KTFEDAQKSI SQHYSKPRVT PVEVMPVFPD FKMWINPCAQ VIFDSDPAPK DTSGAAALEM MSQAMIRGMM DEEGNQFVAY FLPVEETLKK RKR DQEEEMD YAPDDVYDYK IAREYNWNV K NKASKGYEEN YFFIFREGDG VYYNELETRV RLSKRRAKAG VQSGTNALLV VKHRDMNEKE LEAQEARKAQ LENHEPEEEEE EEMETEEKE AGGSDEEQEK GSSSEKEGSE DEHSGSESER EEGDRDEASD KSGSGEDESS EDEARAARDK EEIFGSDADS EDDADS DDED RGQAQGGSDN DSDSGSNGGG QRSRSHRSA SPFPSGSEHS AQEDGSEAAA SDSSEADSDS D
Source	Mammalian cell
Target Names	PAF1
Protein Names	Recommended name: RNA polymerase II-associated factor 1 homolog Short name= hPAF1 Alternative name(s): Pancreatic differentiation protein 2
Expression Region	1-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.