



Recombinant Human Patatin-like phospholipase domain-containing protein 1 (PNPLA1)

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP843295HU |
| Abbreviation | PNPLA1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8N8W4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEEQVFKGDP DTPHSISFSG SGFLSFYQAG AVDALRDLAP RMLETAHRFA GTSAGAVIAA LAICGIEMDE YLRVLNVGVA EVKKSFLGPL SPSCKMVQMM RQFLYRVLPE DSYKVTTGKL HVSLTRLTDG ENVVVSEFTS KEELIEALYC SCFVPVYCGL IPPTYRGVRY IDGGFTGMQP CAFWTDAITI STFSGQQDIC PRDCPAIFHD FRMFNCSFQF SLENIARMTH ALFPPDLVIL HDYYRYGYED AVLYLRRLNA VYLNSSSKRV IFPRVEVYCQ IELALGNECP ERSQPSLRAR QASLEGATQP HKEWVPKGDG RGSHPVVSQ PVQTLFCE SPVSAPVSPL EQPPAQPLAS STPLSLSGMP PVSFPAVHKP PSSTPGSSLP TPPPGLSPLS PQQVQPSGS PARSLHSQAP TSPRPSLGPS TVGAPQTLPR SLSAFPAPQ PVEELGQEQP QAVALLVSSK PKSAVPLVHV KETVSKPYVT ESPAEDSNWV NKVFKKNKQK TSGTRKGFPR HSGSKPSSK VQ |
| Source | Mammalian cell |
| Target Names | PNPLA1 |
| Protein Names | Recommended name: Patatin-like phospholipase domain-containing protein 1 EC= 3.1.1.- |
| Expression Region | 1-532 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.