



Recombinant Human Fibrous sheath-interacting protein 1 (FSIP1)

Product Code	CSB-BP843306HU
Abbreviation	FSIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA03
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDIIKGNLDG ISKPASNSRI RPGSRSSNAS LEVLSTEPGS FKVDTASNLN SGKEDHSESS NTENRRTSND DKQESCSEKI KLAEEGSDDED LDLVQHQIIS ECSDEPKLKE LDSQLQDAIQ KMKKLDKILA KKQRREKEIK KQGLEMRIKL WEEIKSAKYS EAWQSKEEME NTKKFLSLTA VSEETVGP SH EEEDTFSSVF HTQIPPEEYE MQMQKLNKDF TCDVERNESL IKSGKKPFSN TEKIELRGKH NQDFIKRNIE LAKESRNPVV MVDREKKRLV ELLKDLDEKD SGLSSSEG DQ SGWVVPVKGY ELAVTQHQQL AEIDIKLQEL SAASPTISSF SPRLENRNNQ KPDRDGERNM EVTPGEKILR NTKEQRDLHN RLREIDEK LK MMKENVLEST SCLSEEQLKC LLDECILKQK SIIKLSSERK KEDIEDVTPV FPQLSRSIIS KLLNESETKV QKTEVEDADM LESECEASK GYYLTKALTG HNMSEALVTE AENMKCLQFS KDVIISDTKD YFMSKTLGIG R LKRP SFLDD PLYGISV SLS SEDQHLKLSS PENTIADEQE TKDAAEECKE P
Source	Baculovirus
Target Names	FSIP1
Protein Names	Recommended name: Fibrous sheath-interacting protein 1
Expression Region	1-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.