



# Recombinant Human Centrosomal protein POC5 (POC5)

<b>Product Code</b>	CSB-MP843309HU
<b>Abbreviation</b>	POC5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NA72
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSDEEKYSL PVVQNDSSRG SSVSSNLQEE YEELLHYAIV TPNIEPCASQ SSHPKGELVP DVRISTIHI LHSQGNNSEV RETAIEVGKG CDFHISSHSK TDESSPVLSP RKPSHPVMDF FSSHLLADSS SPATNSSHTD AHEILVSDFL VSDENLQKME NVLDLWSSGL KTNIISELSK WRLNFIWHR MEMRKEKEKH AAHLKQLCNQ INELKELQKT FEISIGRKDE VISSLSHAIG KQKEKIELMR TFFHWRIGHV RARQDVYEGK LADQYYQRTL LKKVWKVWRS VVQKQWKDVV ERACQARAE VCIQISNDYE AKVAMLSGAL ENAKAEIQRM QHEKEHFEDS MKKAFMRGVC ALNLEAMTIF QNRNDAGIDS TNNKKEEYGP GVQGKEHSAH LDPSAPP MPL PVTSP LLPSP PAAVGGASAT AVPSAASMTS TRAASASSVH VPVSALGAGS AATAASEEMY VPRVVTSAQQ KAGRTITARI TGRCDFASKN RISSSLAIMG VSPPMSSVVV EKHPVTVQT IPQATAAKYP RTIHPESSSTS ASRSLGTRSA HTQSLTSVHS IKVVD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	POC5
<b>Protein Names</b>	Recommended name: Centrosomal protein POC5 Alternative name(s): Protein of centriole 5 Short name= hPOC5
<b>Expression Region</b>	1-575
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.