



Recombinant Human Uncharacterized protein C1orf172 (C1orf172)

Product Code	CSB-BP843312HU
Abbreviation	C1orf172
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NAX2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRPGHPRPA SGPPRLGPWE RPTCLCLETY DKPPQPPPSR RTRRPDPKDP GHHGPESITF ISGSAEPALE SPTCCLLWRP WVWEWCRAAF CFRRCRDCLQ RCGACVRGCS PCLSTEDSTE GTAEANWAKE HNGVPPSPDR APPSRRDGQR LKSTMGSSFY YPDVKLKGIP VYPYPRATSP APDADSCCKE PLADPPMRH SLPSTFASSP RGSEYYSFH ESDLDPENMG SGSMSSREID VLIFKLTTEL FSVHQIDELA KCTSDTVFLE KTSKISDLIS SITQDYHLDE QDAEGR LVRG IIRISTRKSR ARPQTSEGRS TRAAAPTAAA PDSGHETMVG SGLSQDELTV QISQETTADA IARKLRPYGA PGYPASHDSS FQGTDTDSSG APLLQVYC
Source	Baculovirus
Target Names	KDF1
Protein Names	Recommended name: Uncharacterized protein C1orf172
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.