



Recombinant Human Protein FAM98A (FAM98A)

Product Code	CSB-YP843321HU
Abbreviation	FAM98A
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8NCA5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MECDLMETDI LESLEDLGYK GPLLEDGALS QAVSAGASSP EFTKLCAWLV SELRVLCCKLE ENVQATNSPS EAEEFQLEVS GLLGEMNCPY LSLTSGDVTK RLLIQKNCLL LLTYLISELE AARMLCVNAP PKKAQEGGGS EVFQELKGIC IALGMSKPPA NITMFQFFSG IEKKKCLKETL AKVPPNHVVGK PLLKKPMGPA HWEKIEAINQ AIANEYEVRR KLLIKRLDVT VQSFGWSDRA KSQTEKLAKV YQPKRSVLSP KTTISVAHLL AARQDLSKIL RTSSGSIREK TACAINKVLM GRVPDRGGRP NEIEPPPEM PPWQKRQDGP QQQTGGRGGG RGGYEHSSYG GRGGHEQGGG RGGRRGGYDHG GRGGGRGNKH QGGWTDGGSG GGGGYQDGGY RDSGFQPGGY HGGHSSGGYQ GGGYGGFQTS SSYTGSYQGG GGYQQDNRYQ DGGHHGDRGG GRGGRGRRGG RGGRAGQGGG WGGRRGSQNYH QGGQFEQHFG HGGYQYNHSG FGQGRHYTS
Source	Yeast
Target Names	FAM98A
Protein Names	Recommended name: Protein FAM98A
Expression Region	1-519
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>