



# Recombinant Human F-box only protein 15 (FBXO15)

<b>Product Code</b>	CSB-EP843324HU-B
<b>Abbreviation</b>	FBXO15
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NCQ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MATGRGRILQ QHWLGLQTLR GPSRGGGAAR GRARAFGCRK GPGVKLSAGS AALRCHAGGG QHWESSFSCC SGFLDGMPSE ILLKIFSULD AVSLLCTGCV SRRFYHLAND NFIWIGIYST AFSPARSNWK FNSVEKIAMS MSFLSVQDKE AGYWKKEYIT KQIASVKAAL ADILKPVNPY TGLPVKTKEA LRIFGLGWA IILKEKGGKEY IMEHVDLSIN DTSVTVIWYG KKWPCLASLS TLDLCGMTPV FTDWYKTPK HRLRWHSLIA KYNLSHLTIS TMIGCDRLIR IFCLHPGLLV GVWKKEEELA FVMANLHFHH LVERSTLGSA TIPYELPPHS PFLDDSPYEG LHGYQLHVDL HSGGVFYLCG TFRNLFTKRG NIENGHVKLI VIHLKNNREH LPLIGKVGLS WKTDIFDGC I KSCSMMDVTL LDEHGKPFWC FSSPVCLRSP ATPSDSSSFL GQTYNVDYVD AGRVHVELV WIRETEEYLI VNLVLYLSIA KINHWFGEY
<b>Source</b>	E.coli
<b>Target Names</b>	FBXO15
<b>Protein Names</b>	Recommended name: F-box only protein 15
<b>Expression Region</b>	1-510
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.