



Recombinant Human Adenosine deaminase domain-containing protein 2 (ADAD2)

Product Code	CSB-MP843326HU
Abbreviation	ADAD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NCV1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASASQGADD DGSRRKPRLA ASLQISPQPR PWRPLPAQAAQ SAWGPAPAPA TYRAEGGWPG VSVLRDSGPG AGAGVGELGA ARAWENLGEQ MGKAPRVPVP PAGLSLPLKD PPASQAVSLL TEYAASLGIF LLFREDQPPG PCFPFSVSAE LDGVVCPAGT ANSKTEAKQQ AALSALCYIR SQLENPESPQ TSSRPPLAPL SVENILTHEQ RCAALVSAGF DLLLDERSPY WACKGTVAGV ILEREIPRAR GHVKEIYKLV ALGTGSSCCA GWLEFSGQQ HDCHGLVIAR RALLRFLFRQ LLLATQGGPK GKEQSVLAPQ PGP GPPFTLK PRVFLHLYIS NTPKGAARDI YLPPTSEGLL PHSPMRLQA HVLGQLKPVC YVAPSLCDTH VGCLSASDKL ARWAVLGLGG ALLAHLVSPL YSTSLILADS CHDPPTLSRA IHTRPCLDSV LGPCLPPPYV RTALHLFAGP PVAPSEPTPD TCRGLSLNWS LGDPGIEVVD VATGRVKANA ALGPPSRLCK ASFLRAFHQ ARA V GKPYLL ALKTYEAKA GPYQEARRQL SLLLDQQGLG AWPSKPLV GK FRN
Source	Mammalian cell
Target Names	ADAD2
Protein Names	Recommended name: Adenosine deaminase domain-containing protein 2 Alternative name(s): Testis nuclear RNA-binding protein-like
Expression Region	1-583
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.