



# Recombinant Human Protein FAM46D (FAM46D)

<b>Product Code</b>	CSB-EP843338HU-B
<b>Abbreviation</b>	FAM46D
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NEK8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSEIRFTNLT WDQVITLDQV LDEVIPIHGK GNFPTEVVKP KDIIHVVKDQ LIGQGIIVKD ARLNGSVASY ILASHNGISY KDLDVIFGVE LPGNEEFQVV KDAVLDCLLD FLPKDVKKEK LSPDIMKDAY VQKLVKVCNG HDCWSLISLS NNTGKNLELK FVSSLRRQFE FSVDSFQIVL DPMLDFYSDK NAKLTKESYP VVVAESMYGD FQEAMTHLQH KLICTRKPEE IRGGGLLYC SLLVHGFKPA CMSEIKNLER YMCSRFFIDF PHIEEQQKI ESYLHNHFIG EGMTKYDYLM TLHGVVNEST VCLMSYERRQ ILHLITMMAL KVLGELNILP NTQKVTCFYQ PAPYFAAEAR YPIYVIPEPP PVSFQPYHPL HFRGSNGMS
<b>Source</b>	E.coli
<b>Target Names</b>	FAM46D
<b>Protein Names</b>	Recommended name: Protein FAM46D
<b>Expression Region</b>	1-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	Antibodies against This protein were found only in plasma from cancer patients. While it may be a target for immunotherapy, the function of this gene is unknown.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.