



Recombinant Human Protein FAM46D (FAM46D)

Product Code	CSB-MP843338HU
Abbreviation	FAM46D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEK8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSEIRFTNLT WDQVITLDQV LDEVIPIHGK GNFPTEMEVKP KDIIHVVKDQ LIGQGIIVKD ARLNGSVASY ILASHNGISY KDLDVIFGVE LPGNEEFQVV KDAVLDCLLD FLPKDVKKEK LSPDIMKDAY VQKLVKVCNG HDCWSLISLS NNTGKNLELK FVSSLRRQFE FSVDSFQIVL DPMLDFYSDK NAKLTKESYP VVVAESMYGD FQEAMTHLQH KLICTRKPEE IRGGGLLYC SLLVHGFKPA CMSEIKNLER YMCSRFFIDF PHIEEQQKI ESYLHNHFIG EGMTKYDYLM TLHGVVNEST VCLMSYERRQ ILHLITMMAL KVLGELNILP NTQKVTCFYQ PAPYFAAEAR YPIYVIPEPP PVSFQPYHPL HFRGSNGMS
Source	Mammalian cell
Target Names	FAM46D
Protein Names	Recommended name: Protein FAM46D
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Antibodies against This protein were found only in plasma from cancer patients. While it may be a target for immunotherapy, the function of this gene is unknown.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.