



# Recombinant Human Protein ABHD11 (ABHD11)

<b>Product Code</b>	CSB-EP844049HU-B
<b>Abbreviation</b>	ABHD11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NFV4
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRAGQQLASM LRWTRAWRLP REGLGPHGPS FARVPVAPSS SSGGRGGAEP RPLPLSYRLL DGEAALPAVV FLHGLFGSKT NFNSIAKILA QQTGRRVLTV DARNHGDSPPH SPDMSEYIMS QDLQDLLPQL GLVPCVVVGH SMGGKTAMLL ALQRPELVER LIAVDISPVE STGVSHFATY VAAMRAINIA DELPRSRARK LADEQLSSVI QDMAVRQHLL TNLVEVDGRF VWRVNLDAIT QHLDKILAFP QRQESYLGPT LFLGGNSQF VHPSHHPEIM RLFPPRAQMQT VPNAGHWIHA DRPQDFIAAI RGFLV
<b>Source</b>	E.coli
<b>Target Names</b>	ABHD11
<b>Protein Names</b>	Recommended name: Abhydrolase domain-containing protein 11 EC= 3.-.-. Alternative name(s): Williams-Beuren syndrome chromosomal region 21 protein
<b>Expression Region</b>	1-315
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.