



# Recombinant Human NF-kappa-B inhibitor delta (NFKBID)

<b>Product Code</b>	CSB-MP844092HU
<b>Abbreviation</b>	NFKBID
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NI38
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEAGPWVRSV PPSGPPQFPA VVPGPSLEVA RAHMLALGPQ QLLAQDEEGD TLLHLFAARG LRWAAYAAAE VLQVYRRLDI REHKGKTPLL VAAAANQPLI VEDLLNLGAE PNAADHQGRS VLHVAATYGL PGVLLAVLNS GVQVDLEARD FEGLTPLHTA ILALNVAMRP SDLCPRVLST QARDRLDCVH MLLQMGANHT SQEIKSNKTV LHLAVQAANP TLVQLLLELP RGDRLTFVNM KAHGNTALHM AAALPPGPAQ EAIVRHLLAA GADPTLRNLE NEQPVHLLRP GPGPEGLRQL LKRSRVAPPG LSS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NFKBID
<b>Protein Names</b>	Recommended name: NF-kappa-B inhibitor delta Alternative name(s): I-kappa-B-delta Short name= IκB-delta Short name= IκappaBdelta IκappaBNS T-cell activation NFKB-like protein TA-NFKBH
<b>Expression Region</b>	1-313
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.