



# Recombinant Mouse Abhydrolase domain-containing protein 8 (Abhd8)

<b>Product Code</b>	CSB-MP844428MO
<b>Abbreviation</b>	Abhd8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8R0P8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSEHTMLTG VTDGFFCCLL GTPPNAVRPL ESVESSDGYT FVEVKPGRVL RVKHAGPAPI PTPPPPPED DPGVKTGLVR CQRRITVYRN GRLVVENLGR APRADLQGRS GSGDPPAALE VELAEPAGGD TRANPGSGRR RRRRRPKRTI HIDCEQRITS CKGAQADVVL FFIHGVGGSL AIWKEQLDFF VRLGYEVVAP DLAGHGASSA PQVAAAAYTFY ALAEDMRAIF TRYAKKRNVL IGHSYGVSFC TFLAHEYPLD VHKVIMINGG GPTALEPSLC SIFNMPTCVL HCLSPCLAWS FLKAGFARQG AKEKQLLKEG NAFNVSSFVL RAMMSGQYWP EGDEVYHAEL TVPVLLVHGM HDKFVPVEED QRMAEILLLA FLKLIEEGSH MVMLECPETV NTLLHEFLLW EPEPEAEPKL EPKPKPQLLQ PEPAPGEEK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Abhd8
<b>Protein Names</b>	Recommended name: Abhydrolase domain-containing protein 8 EC= 3.-.-.
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. The predicted protein contains alpha/beta hydrolase fold and secretory lipase domains.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.