



Recombinant Mouse Amyloid beta A4 precursor protein-binding family B member 1-interacting protein (Apbb1ip)

Product Code	CSB-BP844466MO
Abbreviation	Apbb1ip
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R5A3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MGESNEDIDQ MFSTLLGEMD LLTQSLGVDT LPPDPNPPR EEFNYTVGFK DLNESLNALE DQDLDALMAD LVADISEAEQ RTIQAQKESS QNQDRFALLR ASDGGGTASG GYGASAAID VSHHEEALPP PPVEPMLDLL PPPPPPPPE LLSKEEEAK AKADKIKLAL EKLKEAKVKK LVVKVHMDDS STKSLMVDER QLARDVLDNL FEKTHCDCNV DWCLYEIYPE LQIERVFEDH ENVVEVLSDW TRDTENKVLV LEKEERYAVF KNPQNFYLDN KGKKNKETS EKMNANKNEY LLEESFCGTS IIVPELEGAL YLKEDGKKSX KRRYFLLRAS GIYYVPGKGT KTSRDLACFI QFENVNIYYG IQCKMKYKAP TDHCFVLKHP QIQKESQYIK YLCCDDARTL SQWVMGIRIA KYGKTLYDNY QRAVARAGLA SRWTNLGTVG TPMPAQSTV SSSLKGTGTSQ PNGQMPQAIP SAGPPLQEAQ TQIETTKDEK QGLGNHSPGA TRENHRPKSS LPPPPPPVRR SSDTCGSPAL PSKVKGPGTC TFPHPPENFL PPPPPPPPEE DNSGLLPPPP PPPYLEEPPD FVPPPPPPAA VEDSALPPPP PPPPCLSQEI TKSSPLPPPK PLVPPKRQEN QGLPGAPGNS EQDFMSDLMK ALQKKRGNIP
Source	Baculovirus
Target Names	Apbb1ip
Protein Names	Recommended name: Amyloid beta A4 precursor protein-binding family B member 1-interacting protein Alternative name(s): APBB1-interacting protein 1 Proline-rich EVH1 ligand 1 Short name= PREL-1 Proline-rich protein 48
Expression Region	1-670
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.