



# Recombinant Mouse Amyloid beta A4 precursor protein-binding family B member 1-interacting protein (Apbb1ip)

<b>Product Code</b>	CSB-YP844466MO
<b>Abbreviation</b>	Apbb1ip
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8R5A3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGESNEDIDQ MFSTLLGEMD LLTQSLGVDT LPPDPNPPR EEFNYTVGFK DLNESLNALE DQDLDALMAD LVADISEAEQ RTIQAQKESS QNQDRFALLR ASDQQGTASG GYGASAAID VSHHEEALPP PPVEPMLDLL PPPPPPPPE LLSKEEEEEAK AKADKIKLAL EKLKEAKVKK LVVKVHMDDS STKSLMVDER QLARDVLDNL FEKTHCDCNV DWCLYEIYPE LQIERVFEDH ENVVEVLSDW TRDTENKVLV LEKEERYAVF KNPQNFYLDN KGKKNKETS EKMNANKNEY LLEESFCGTS IIVPELEGAL YLKEDGKKS W KRRYFLLRAS GIYYVPGKGT KTSRDLACFI QFENVNIYYG IQCKMKYKAP TDHCFVLKHP QIQKESQYIK YLCCDDARTL SQWVMGIRIA KYGKTLYDNY QRAVARAGLA SRWTNLGTVG TPMPAQSTV SSSLKGTGTSQ PNGQMPQAIP SAGPPLQEAQ TQIETTKDEK QGLGNHSPGA TRENHRPKSS LPPPPPPVRR SSDTCGSPAL PSKVKGPGTC TFPHPPENFL PPPPPPPPEE DNSGLLPPPP PPPYLEEPPD FVPPPPPPAA VEDSALPPPP PPPPCLSQEI TKSSPLPPK PLVPPKRQEN QGLPGAPGNS EQDFMSDLMK ALQKKRGNIP
<b>Source</b>	Yeast
<b>Target Names</b>	Apbb1ip
<b>Protein Names</b>	Recommended name: Amyloid beta A4 precursor protein-binding family B member 1-interacting protein Alternative name(s): APBB1-interacting protein 1 Proline-rich EVH1 ligand 1 Short name= PREL-1 Proline-rich protein 48
<b>Expression Region</b>	1-670
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.