



# Recombinant Goat Chitinase-3-like protein 1 (CHI3L1)

<b>Product Code</b>	CSB-BP844620GO
<b>Abbreviation</b>	CHI3L1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8SPQ0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Capra hircus (Goat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	YKLCYYTS WSQYREGDGS CFPDAIDPFL CTHIIYSFAN ISNNEIDTWE WNDVTLYDTL NTLKRNPKL KTLLSVGGWN FGPERFSKIA SKTQSRRTFI KSVPPFLRTH GFDGLDLAWL YPGRDRKRHL TGLVKEMKAE FAREAQAGTE RLLSAAVSA GKIAIDRGYD IAQISRHLDF ISLLTYDFHG AWRQTVGHHS PLFRGQEDAS SDRFSNADYA VSYMLRLGAP ANKLVMGIPT FGRSFTLASS KTDVGAPISG PGIPGRFTKE KGILAYYEIC DFLHGATTHR FRDQQVPYAT KGNQWVAYDD QESVKNKARY LKNRQLAGAM VWALDLDDFR GTFCGQNLTF PLTSAVKDVL AEV
<b>Source</b>	Baculovirus
<b>Target Names</b>	CHI3L1
<b>Protein Names</b>	Recommended name: Chitinase-3-like protein 1 Alternative name(s): BP40 Mammary gland protein 40 Short name= MGP-40
<b>Expression Region</b>	22-383
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Chitinases catalyze the hydrolysis of chitin, which is an abundant glycopolymer found in insect exoskeletons and fungal cell walls. The glycoside hydrolase 18 family of chitinases includes eight human family members. This gene encodes a glycoprotein member of the glycosyl hydrolase 18 family. The protein lacks chitinase activity and is secreted by activated macrophages, chondrocytes, neutrophils and synovial cells. The protein is thought to play a role in the process of inflammation and tissue remodeling.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.