



Recombinant Human Smad nuclear-interacting protein 1 (SNIP1)

Product Code	CSB-MP844702HU
Abbreviation	SNIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAD8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKAVKSERER GSRRRHRDGD VVLPAGVVVK QERLSPEVAP PAHRRPDHSG GSPSPPTSEP ARSGHRGNRA RGVSRSPPKK KNKASGRRSK SPRSKRNRSP HHSTVKVKQE REDHPRRGRE DRQHREPSEQ EHRRARNSDR DRHRGHSHQR RTSNERPGSG QGQGRDRDTQ NLQAQEEERE FYNARRREHR QRNDVGGGGS ESQELVPRPG GNNKEKEVPA KEKPSFELSG ALLEDTNTFR GVVIKYSEPP EARIPKKRWR LYPFKNDEVL PVMYIHRQSA YLLGRHRIA DIPIDHPSCS KQHAVFYRL VEYTRADGTV GRRVKPYIID LGSGNGTFLN NKRIEQRYY ELKEKDVLKF GFSSREYVLL HESSDTSEID RKDDEDEEEE EEVSDS
Source	Mammalian cell
Target Names	SNIP1
Protein Names	Recommended name: Smad nuclear-interacting protein 1 Alternative name(s): FHA domain-containing protein SNIP1
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.