



Recombinant Human GA-binding protein subunit beta-2 (GABPB2)

Product Code	CSB-MP844704HU
Abbreviation	GABPB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAK5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLVDLGKRL LEAARKGQDD EVRTL MANGA PFTTDWLGTS PLHLAAQYGH YSTAEVLLRA GVS RDARTKV DRTPLHMAAA DGH AHIVELL VRNGADVNAK DMLKMTALHW ATERHHRDVV ELLIKYGADV HAFSKFDKSA FDIALEKNN A EILVILQEAM QNQVNVNPER ANPVTDPVSM AAPFIFTSGE VVNLASLISS TNTKTTSGDP HASTVQFSNS TTSVLATLAA LAEASVPLSN SHRATANTEE IIEGNSVDSS IQQVMGSGGQ RVITIVTDGV PLGNIQTSIP TGGIGQPFIV TVQDGQQVLT VPAGKVAEET VIKEEEEELK PLTKKPRIGE KTNSVEESKE GNERELLQQQ LQEANRRAQE YRHQLLKEEQ EAEQYRLKLE AIARQQPNGV DFTMVEEVAE VDAVVVTEGE LEERETKVTG SAGTTEPHTR VSMATVSS
Source	Mammalian cell
Target Names	GABPB2
Protein Names	Recommended name: GA-binding protein subunit beta-2 Short name= GABP subunit beta-2 Short name= GABPB-2
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.